

## CD45RO

### IN VITRO DIAGNOSTIC DATASHEET

**INTENDED USE :** IN VITRO DIAGNOSTIC USE

This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.

**DESCRIPTION :**

CD45RO is expressed on most thymocytes, monocytes, macrophages and granulocytes. It is not expressed on normal B lymphocytes.

**CATALOG NO :** PL112

PL112-R7 7 ML RTU 70 TEST

PL112-R1 1 ML RTU 10 TEST

**STAINING PATTERN :** Cell Membrane

PL112-1 1 ML 1/200 2000 TEST

PL112-0,1 0,1 ML 1/200 200 TEST

**POSITIVE CONTROL :** Tonsil or lymph node

**VOLUME :** 7 ml Ready to Use ( 7 ml of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing stabilizing protein and 0.015mol/L sodium azide. )

**HOST :** Mouse

**CLONE :** UCHL-1

**ANTIBODY CONCENTRATION :** 100ug / ml

**SPECIES REACTIVITY :** Human, Others not tested.

**EPITOPE :** Not determined

**MICROBIOLOGICAL STATE :** This product is not sterile.

**PRETREATMENT :** No special pretreatment is required for immuno-histochemistry of formalin-fixed tissues.

**PRIMARY ANTIBODY INCUBATION TIME :** 30 minutes at Room Temperature

**STAINING TIPS :** If the staining is too light, use lower dilution or longer time. If the staining is too strong, check pretreatment, use higher dilution or shorter time.

**STORAGE AND STABILITY :** This product contains sodium azide and is stable for 24 months when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent is not stored as recommended, performance must be validated by the user.

**TROUBLESHOOTING :** Please contact Patolab Technical Support by e-mail ( patolab@patolab.com.tr ).