

MAYER HEMATOXYLIN

INTENDED USE:

For In Vitro Diagnostic Use

<u>Catalog No</u>	<u>Volume</u>
PL-125-MH	125 ml
PL-1000-MH	1000 ml

DESCRIPTION

Mayer's Hematoxylin is a stain with histological applications. This formula is based on the original formulation developed by Mayer and requires slightly longer application periods.

WARNINGS & PRECAUTIONS

Refer to MSDS.

STORAGE & SHELF LIFE

Store at room temperature. Each component is stable for 24 months.

MICROBIOLOGICAL STATE

Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED

N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE

1. Rehydrate tissue if needed.
2. Stain nuclei 1-2 minutes with hematoxylin.
3. Place in Bluing Reagent until stain is blue.
4. Rinse in DI water.
5. Coverslip.

Troubleshooting

Please contact PatoLab Technical Support by e-mail (patolab@patolab.com.tr).