

## Calcitonin

IN VITRO DIAGNOSTIC DATASHEET

**INTENDED USE :** IN VITRO DIAGNOSTIC USE

This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy. It is Rabbit Monoclonal.

**DESCRIPTION :** Calcitonin is a 32 amino acid peptide which can be demonstrated in C cells of the normal and hyperplastic thyroid. Staining for calcitonin may be used for the identification of a spectrum of C cell proliferative abnormalities ranging from C cell hyperplasia to invasive tumors. Staining for calcitonin in medullary carcinoma of the thyroid produces a fine granular pattern in the cytoplasm. Amyloid deposits within the tumor may also exhibit varying degrees of calcitonin activity.

<b>CATALOG NO :</b>	PL9117	PL9117-R7	7 ML RTU 70 TEST
		PL9117-R1	1 ML RTU 10 TEST
<b>STAINING PATTERN :</b>	Cytoplasmic	PL9117-1	1 ML 1/100 1000 TEST
		PL9117-0,1	0,1 ML 1/100 100 TEST

**POSITIVE CONTROL :** Thyroid or medullary CA

**VOLUME :** 7 ml Ready to Use ( 7 ml of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing stabilizing protein and 0.015mol/L sodium azide. )

**HOST :** Rabbit

**CLONE :** SP17

**ANTIBODY CONCENTRATION :** Not known

**SPECIES REACTIVITY :** Human, Monkey, Mouse, Rat, Dog, Horse, Sheep and Pig. Others not tested.

**EPITOPE :** 1-32 aa

**MICROBIOLOGICAL STATE :** This product is not sterile.

**PRETREATMENT :** No special pretreatment is required for immuno-histochemistry of formalin-fixed tissues.

**PRIMARY ANTIBODY INCUBATION TIME :** 30 minutes at Room Temperature

**STAINING TIPS :** If the staining is too light, use lower dilution or longer time. If the staining is too strong, check pretreatment, use higher dilution or shorter time.

**STORAGE AND STABILITY :** This product contains sodium azide and is stable for 24 months when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent is not stored as recommended, performance must be validated by the user.

**TROUBLESHOOTING :** Please contact Patolab Technical Support by e-mail ( patolab@patolab.com.tr ).