

CD8

IN VITRO DIAGNOSTIC DATASHEET

INTENDED USE: IN VITRO DIAGNOSTIC USE

This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy. it is Rabbit Monoclonal.

DESCRIPTION: CD8 molecule consists of two chains, termed alpha and Beta chain, which are expressed as a disulphide-linked alpha/Beta heterdimer or as an alpha/alpha homodimer on T cell subset, thymocytes and NK cells. The majority of CD8+ T cells express CD8 as alpha/Beta heterdimer. CD8 functions as a coreceptor in concert with TCR for binding the MHC class I/peptide complex. The HIV-2 envelope glycoprotein binds CD8 alpha chain (but not Beta chain).

 CATALOG NO:
 PL9116
 PL9116-R7
 7 ML RTU 70 TEST

 PL9116-R1
 1 ML RTU 10 TEST

 STAINING PATTERN:
 Cell Membrane
 PL9116-1
 1 ML 1/500 5000 TEST

 PL9116-0,1
 0,1 ML 1/500 500 TEST

POSITIVE CONTROL: Tonsil

VOLUME: 7 ml Ready to Use (7 ml of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing

stabilizing protein and 0.015mol/L sodium azide.)

HOST: Rabbit

CLONE: SP16

ANTIBODY CONCENTRATION: Not known

SPECIES REACTIVITY: Human. Others not tested.

EPITOPE: C-terminus

MICROBIOLOGICAL STATE: This product is not sterile.

PRETREATMENT: Staining of formalin-fixed tissue sections requires treating the tissue sections in boiling 10mM Citrate, pH 6.0, for 10-20 min followed by cooling at RT for 20 min.

PRIMARY ANTIBODY INCUBATION TIME: 30 minutes at Room Temperature

STAINING TIPS: If the staining is too light, use lower dilution or longer time. If the staining is

too strong, check pretreatment, use higher dilution or shorter time.

STORAGE AND STABILITY: This product contains sodium azide and is stable for 24 months when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent

is not stored as recommended, performance must be validated by the user.

TROUBLESHOOTING: Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).



