

Harris Hematoxylin

INTENDED USE: For In Vitro Diagnostic Use

<u>Catalog No</u>	<u>Volume</u>
PL-500-HAR	500 ml
PL-1000-HAR	1000 ml

DESCRIPTION

Nuclear stain for Papanicolaou method.

Product for the preparation of: gynecological specimens, urine cytology, fine needle specimens, sputum and bronchial washings, to be examined by optical microscopy.

STORAGE & SHELF LIFE : Store at room temperature. Stable for minimum 3 years.

MICROBIOLOGICAL STATE : Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED : N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE

- 1) Ethanol 95°, 2 minutes
- 2) Distilled water, 2 minutes
- 3) Harris Hematoxylin, 1 minute
- 4) Tap water, 5 minutes
- 5) Ethanol 95°, 15 seconds
- 6) OG 6, 2 minutes
- 7) Ethanol 95°, 15 seconds (twice)
- 8) EA 65, 5 minutes
- 9) Ethanol 95°, 15 seconds
- 10) Absolute Ethanol, 30 seconds (twice)
- 11) Xylene, 2 minutes (twice)

RESULTS

Nuclei : Blue-purple
Cyanophil cytoplasm : Blue-green
Eosinophil cytoplasm : Pink
Keratinized cytoplasm : From pink to orange

Troubleshooting

Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).