

IN VITRO DIAGNOSTIC DATASHEET

CD4

INTENDED USE : IN VITRO DIAGNOSTIC USE

PatoLab Biyomedikal Cihazlar

patolab@patolab.com.tr

This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.

DESCRIPTION :

CD4, a single chain transmembrane glycoprotein, is found on a T cell subset (helper/inducer) r 45% of peripheral blood lymphocytes. It is also present on 80% of thymocytes and at a lower level on monocytes. It is involved in recognition of antigen presented along with MHC class II by APCs. It serves as receptor for HIV.

CATALOG NO : PL1528	PL1528-R7	7 ML RTU 70 TEST
	PL1528-R1	1 ML RTU 10 TEST
STAINING PATTERN : Cell Membrane	PL1528-1	1 ML 1/100 1000 TEST
	PL1528-0,1	0,1 ML 1/100 100 TEST
POSITIVE CONTROL : Tonsil		
VOLUME : 7 ml Ready to Use (7 ml of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing		
	g protein and 0.015mol/L sodium	azide.)
HOST : Mouse		
CLONE : 4B12		
ANTIBODY CONCENTRATION : Not known		
SPECIES REACTIVITY : Human. Others not-tested		
EPITOPE : Not determined		
EFITOPE . Not determined		
MICROBIOLOGICAL STATE : This product is not sterile.		
PRETREATMENT : Staining of formalin-fixed tissue sections requires treating the tissue sections in boiling		
1mM EDTA, pH 8.0, for 10-20 min f	ollowed by cooling at RT for	20 min.
REMARY ANTIDORY INCURATION TIME		
PRIMARY ANTIBODY INCUBATION TIME	: 60 minutes at Room Temper	ature
STAINING TIPS : If the staining is too light, use lower dilution or longer time. If the staining is		
too strong, check pretreatment, use higher dilution or shorter time.		
STORAGE AND STABILITY : This product	contains sodium azide and is stab	le for 24 months when stored
	fter expiration date indicated on la	
is not stored as recom	mended, performance must be va	lidated by the user.
	Datalah Tashaisal Osma art burya	
TROUBLESHOOTING : Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).		

www.patolab.com.tr

Sarf

Mal

a n

IVD