Pat Lab[®]

ANTI-POLYVALENT HRP KIT

INTENDED USE:

For In Vitro Diagnostic Use

CATALOG NO	VOLUME		
PL-060-HL PL-125-HL PL-1000-HL	60 ml (600 slides) 125 ml (1250 slides) 1000 ml (10000 slides)		
SPECIFICITY :	Anti-Mouse IgG , Anti-Rabbit IgG		
ENZYME :	Peroxidase		

REAGENTS :

Quantity	Component	PL-060-HL	PL-125-HL	PL-1000-HL
1	Protein Block	PL-060-PBQ	PL-125-PBQ	PL-1000-PBQ
1	Biotinylated Goat Anti-Polyvalent	PL-060-BN	PL-125-BN	PL-1000-BN
1	Streptavidin Peroxidase	PL-060-HR	PL-125-HR	PL-1000-HR

DESCRIPTION

The reagents in this kit constitute a labeled streptavidin-biotin immunoenzymatic antigen detection system. This technique involves the sequential incubation of the specimen with an unconjugated primary antibody specific to the target antigen, a biotinylated secondary antibody that reacts with the primary antibody, enzyme-labeled streptavidin, and substrate chromogen.

PRINCIPLE OF THE PROCEDURE

This HRP detection system detects a specific antibody bound to an antigen in tissue sections. The specific antibody is located by a biotin-conjugated secondary antibody. This step is followed by the addition of a streptavidin-enzyme conjugate that binds to the biotin present on the secondary antibody. The specific antibody, secondary antibody, and streptavidin-enzyme complex is then visualized with an appropriate substrate/chromogen.

WARNINGS & PRECAUTIONS

Refer to MSDS.

STORAGE & SHELF LIFE

Store at 2-8°C. Product is stable for 24 months from the date of manufacture.

MICROBIOLOGICAL STATE

Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED

Primary antibody.

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PatoLab Biyomedikal Cihazlar ve Sarf Malzemeler San. ve Tic. Ltd. Şti. patolab@patolab.com.tr www.patolab.com.tr



Pat Lab[®]

PROCEDURE

STAINING PROTOCOL (kit components in bold):

- 1. Deparaffinize and rehydrate tissue section.
- 2. To reduce nonspecific background staining due to endogenous peroxidase, incubate slide in hydrogen peroxide for 10-15 minutes.
- 3. Wash 3 times in buffer.
- 4. If required, incubate tissue in digestive enzyme (or appropriate pretreatment).
- 5. Wash 3 times in buffer.
- 6. (Optional) Apply **Protein Block** and incubate for 5 minutes at room temperature to block nonspecific background staining.
 - NOTE: Do not exceed 10 minutes or there may be a reduction in desired stain.
- 7. Rinse (Optional).
- 8. Apply primary antibody and incubate according to manufacturer's protocol.
- 9. Wash 4 times in buffer.
- 10. Apply Biotinylated Goat Anti-Polyvalent and incubate for 10 minutes at room temperature.
- 11. Wash 3 times in buffer.
- 12. Apply Streptavidin Peroxidase and incubate for 10 minutes at room temperature.
- 13. Rinse 3 times in buffer.
- 14. Incubate with peroxidase-compatible chromogen of choice according to manufacturer's recommendations.
- 15. Counterstain and coverslip.

The specificity and sensitivity of antigen detection is dependent on the specific primary antibody used.

Troubleshooting

Please contact PatoLab Technical Support by e-mail (patolab@patolab.com.tr).



IVD