

Rapid Decalcifier

INTENDED USE : For In Vitro Diagnostic Use

<u>Catalog No</u>	<u>Volume</u>
PL-050-RDC	50 ml
PL-1000-RDC	1000 ml

DESCRIPTION

Rapid Decalcifier, based on hydrochloric acid. It acts on all mineralized tissues: compact bone, concretions or calcareous deposits in fleshy tissue. A saline corrector has been added to check tissue swelling caused by acid substances.

STORAGE & SHELF LIFE

Store at room temperature. Stable for minimum 3 years.

MICROBIOLOGICAL STATE

Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED

N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE

- 1) Volume ratio specimen/decalcifier 1 : 100
- 2) Procedure time 4 - 8 hours - < 5 mm thickness.
- 3) Post decalcifying procedure: reconditioning of the sample. Rinse with slowly running water for 1 hour or apply 3 changes of PBS pH 7.4, 20 minutes each.

Troubleshooting

Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).