

## S100 Protein

## IN VITRO DIAGNOSTIC DATASHEET

## INTENDED USE : IN VITRO DIAGNOSTIC USE

This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.

## **DESCRIPTION**:

S100 belongs to the family of calcium binding proteins such as calmodulin and troponin C. S100A is composed of an alpha and beta chain whereas S100B is composed of two beta chains. S100 protein is also expressed in the antigen presenting cells such as the Langerhans cells in skin and interdigitating reticulum cells in the paracortex of lymph nodes.

| CATALOG NO : PL296  | PL296-R7            | 7 ML RTU 70 TEST                         |
|---|---------------------|--|
| STAINING PATTERN : Cytoplasmic  | PL296-R1<br>PL296-1 | 1 ML RTU 10 TEST<br>1 ML 1/100 1000 TEST |
|   | PL296-0,1           | 0,1 ML 1/100 100 TEST                    |
| <b>POSITIVE CONTROL</b> : Melanoma or Schwannoma  |                     |  |
| <b>VOLUME:</b> 7 ml Ready to Use ( 7 ml of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing stabilizing protein and 0.015mol/L sodium azide. )  |                     |  |
| HOST : Mouse  |                     |  |
| <b>CLONE :</b> 4C4.9  |                     |  |
| ANTIBODY CONCENTRATION : 200ug / ml   |                     |  |
| SPECIES REACTIVITY: Human, Cow, Rat, and Mouse. Others-not known  |                     |  |
| EPITOPE : Belirlenmedi.   |                     |  |
| MICROBIOLOGICAL STATE : This product is not sterile.  |                     |  |
| <b>PRETREATMENT</b> : No special pretreatment is required for immuno-histochemistry of formalin-fixed tissues.  |                     |  |
| <b>PRIMARY ANTIBODY INCUBATION TIME :</b> 30 minutes at Room Temperature  |                     |  |
| <b>STAINING TIPS :</b> If the staining is too light, use lower dilution or longer time. If the staining is too strong, check pretreatment, use higher dilution or shorter time.   |                     |  |
| <b>STORAGE AND STABILITY :</b> This product contains sodium azide and is stable for 24 months when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent is not stored as recommended, performance must be validated by the user. |                     |  |
| TROUBLESHOOTING: Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).   |                     |  |

PatoLab Biyomedikal Cihazlar ve Sarf Malzemeler San. ve Tic. Ltd. Şti.patolab@patolab.com.trwww.patolab.com.tr

