

## FITC MOUNTING MEDIA

### INTENDED USE:

For In Vitro Diagnostic Use

<u>Catalog No</u>	<u>Volume</u>
PL-060-FM	60 ml
PL-125-FM	125 ml

### DESCRIPTION

Immobilizes coverglass over section, enhances FITC intensity and reduces fading to a significant degree

### WARNINGS & PRECAUTIONS

Reagent should be at room temperature before opening vial for use. **DO NOT FREEZE.**  
**NO MSDS IS REQUIRED FOR THIS PRODUCT.** The only chemical of a toxic nature in this product is sodium azide, which is present at a concentration below 0.1%.

### STORAGE & SHELF LIFE

Store at room temperature, stable for 24 months.

### MICROBIOLOGICAL STATE

Product(s) not sterile.

### MATERIALS REQUIRED BUT NOT PROVIDED

N/A

### SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

### PROCEDURE

#### **FOR USE IN FLUORESCENCE MICROSCOPY**

After completion of immunostaining, remove slide from water, apply sufficient drops of FITC Mounting Media upon section and mount coverglass. Coverglass should be immobilized within 30 minutes at room temperature. To avoid air bubbles within mounting media, DO NOT SHAKE VIAL. If air bubbles are evident during coverslipping, to avoid fingerprint, push the coverslip with your nail. These directions are recommendations; however, they may be altered according to procedures established in your laboratory.

### Troubleshooting

Please contact PatoLab Technical Support by e-mail ([patolab@patolab.com.tr](mailto:patolab@patolab.com.tr)).