

Papanicolaou OG6

INTENDED USE : For In Vitro Diagnostic Use

Catalog No	Volume
PL-500-OG6	500 ml
PL-1000-OG6	1000 ml

DESCRIPTION

Product for the preparation of: gynecological specimens, urine cytology, fine needle specimens, sputum and bronchial washings, to be examined by optical microscopy

STORAGE & SHELF LIFE

Store at room temperature. Stable for minimum 3 years.

MICROBIOLOGICAL STATE : Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED : N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE

- 1) Ethanol 95° 2 minutes
- 2) Distilled water 2 minutes
- 3) Harris Hematoxylin 1 minute
- 4) Tap water 5 minutes
- 5) Ethanol 95° 15 seconds
- 6) OG 6 2 minutes
- 7) Ethanol 95° 15 seconds (twice)
- 8) EA 50 5 minutes
- 9) Ethanol 95° 15 seconds
- 10) Absolute Ethanol 30 seconds (twice)
- 11) Xylene 2 minutes (twice)

RESULTS

Nuclei	Blue-Purple
Cyanophil cytoplasm	Blue-green
Eosinophil cytoplasm	Pink
Keratinized cytoplasm	From Pink to orange

Troubleshooting

Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).