

CD31 / PECAM-1

IN VITRO DIAGNOSTIC DATASHEET

INTENDED USE : IN VITRO DIAGNOSTIC USE

This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.

DESCRIPTION:

CD31 is a glycoprotein expressed on endothelial cells and in platelets. It is known to be involved in cell signaling and cell adhesion.

| CATALOG NO : PL353 | PL353-R7 | 7 ML RTU 70 TEST |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|-----------------------|
| | PL353-R1 | 1 ML RTU 10 TEST |
| STAINING PATTERN : Cell Membrane | PL353-1 | 1 ML 1/100 1000 TEST |
| | PL353-0,1 | 0,1 ML 1/100 100 TEST |
| POSITIVE CONTROL : Placenta | | |
| VOLUME : 7 ml Ready to Use (7 ml of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing stabilizing protein and 0.015mol/L sodium azide.) | | |
| HOST: Mouse | | |
| CLONE : JC/70A | | |
| ANTIBODY CONCENTRATION : Not known | | |
| SPECIES REACTIVITY : Human. Does not react with rat and pig. Others-not tested. | | |
| EPITOPE : Not determined | | |
| MICROBIOLOGICAL STATE : This product is not sterile. | | |
| PRETREATMENT : Staining of formalin-fixed tissue sections requires treating the tissue sections in boiling 1mM EDTA, pH 8.0, for 10-20 min followed by cooling at RT for 20 min. | | |
| PRIMARY ANTIBODY INCUBATION TIME : 30 minutes at Room Temperature | | |
| STAINING TIPS : If the staining is too light, use lower dilution or longer time. If the staining is too strong, check pretreatment, use higher dilution or shorter time. | | |
| STORAGE AND STABILITY : This product contains sodium azide and is stable for 24 months when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent is not stored as recommended, performance must be validated by the user. | | |
| TROUBLESHOOTING : Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr). | | |

PatoLab Biyomedikal Cihazlar ve Sarf Malzemeler San. ve Tic. Ltd. Şti.patolab@patolab.com.trwww.patolab.com.tr

