

## DOG-1

## IN VITRO DIAGNOSTIC DATASHEET

INTENDED USE: IN VITRO DIAGNOSTIC USE

This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.

## **DESCRIPTION:**

Dog1 is a cell surface protein of unknown function selectively expressed in gastrointestinal stromal tumors (GIST). Among GIST cases with KIT mutations the Dog1 antibody identified 11% more cases than c-Kit antibody. Dog1 identifies vast majority of both c-Kit negative and PDGFRA mutated GIST cases that may still benefit from imatinib mesylate (Gleevac), an inhibitor of the kit tyrosin kinase.

**CATALOG NO:** PL1933 PL1933-R7 7 ML RTU 70 TEST

PL1933-R1 1 ML RTU 10 TEST

**STAINING PATTERN:** Cytoplasmic, Membranous PL1933-1 1 ML 1/200 2000 TEST

> PL1933-0,1 0,1 ML 1/200 200 TEST

**POSITIVE CONTROL: Gastrointestinal Stromal Tumors** 

**VOLUME:** 7 ml Ready to Use (7 ml of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing

stabilizing protein and 0.015mol/L sodium azide.)

HOST: Mouse

**CLONE:** 1.1

**ANTIBODY CONCENTRATION:** Not known

SPECIES REACTIVITY: Human. Others not-tested

**EPITOPE:** Not determined

**MICROBIOLOGICAL STATE:** This product is not sterile.

PRETREATMENT: Staining of formalin-fixed tissue sections requires treating the tissue sections in boiling 10mM citrate buffer, pH 6.0, for 10-20 minutes followed by cooling at room temperature for 20 min.

PRIMARY ANTIBODY INCUBATION TIME: 60 minutes at Room Temperature

**STAINING TIPS:** If the staining is too light, use lower dilution or longer time. If the staining is

too strong, check pretreatment, use higher dilution or shorter time.

STORAGE AND STABILITY: This product contains sodium azide and is stable for 24 months when stored

at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent

is not stored as recommended, performance must be validated by the user.

TROUBLESHOOTING: Please contact Patolab Technical Support by e-mail ( patolab@patolab.com.tr ).



