

c-erbB-2 / HER-2

IN VITRO DIAGNOSTIC DATASHEET

INTENDED USE : IN VITRO DIAGNOSTIC USE

This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.

DESCRIPTION :

c-erbB-2 is a receptor tyrosine kinase of the c-erbB family. Approximately 25% of primary breast and ovarian tumors were found to overexpress the protein.

CATALOG NO : PL730 STAINING PATTERN : Cell membrane	PL730-R7 PL730-R1 PL730-1 PL730-0,1	7 ML RTU 70 TEST 1 ML RTU 10 TEST 1 ML 1/400 4000 TEST 0,1 ML 1/400 400 TEST	
POSITIVE CONTROL : LS174T cells, or breast CA	FL730-0,1	0,1 ME 1/400 400 TEST	
VOLUME: 7 ml Ready to Use (7 ml of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing stabilizing protein and 0.015mol/L sodium azide.) HOST: Mouse			
CLONE : e2-4001 + 3B5			
ANTIBODY CONCENTRATION : 200ug/ml			
SPECIES REACTIVITY : Human. Others not tested.			
EPITOPE : Cytoplasmic			
MICROBIOLOGICAL STATE : This product is not sterile.			
PRETREATMENT : Staining of formalin-fixed tissue sections requires treating the tissue sections in boiling 10mM citrate buffer, pH 6.0, for 10-20 minutes followed by cooling at room temperature for 20 min.			
PRIMARY ANTIBODY INCUBATION TIME : 30 minutes at Room Temperature			
STAINING TIPS : If the staining is too light, use lower di too strong, check pretreatment, use h	ver dilution or longer time. If the staining is use higher dilution or shorter time.		
STORAGE AND STABILITY : This product contains sodium azide and is stable for 24 months when stored at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent is not stored as recommended, performance must be validated by the user.			
TROUBLESHOOTING · Dease contact Patolab Technical Support by e-mail (natolab@natolab.com tr)			

TROUBLESHOOTING: Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).

