

## MLH<sub>1</sub>

## IN VITRO DIAGNOSTIC DATASHEET

INTENDED USE: IN VITRO DIAGNOSTIC USE

This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy. It is Rabbit Monoclonal.

**DESCRIPTION:** Hereditary non-polyposis colorectal cancer (HNPCC) or Lynch syndrome is an autosomal dominant disease primarily associated with germline mutations in DNA mismatch repair genes. HNPCC accounts for 2-3% of all colorectal cancer. Approximately 70-80% of HNPCC cases result from defects in MLH1 and MSH2. Mutations in MLH1 affect the MLH1-PMS2 complex formation and may cause defects in the DNA error repair mechanism.

 CATALOG NO:
 PL1941
 PL1941-R7
 7 ML RTU 70 TEST

 PL1941-R1
 1 ML RTU 10 TEST

 STAINING PATTERN:
 Nuclear
 PL1941-1
 1 ML 1/100 1000 TEST

 PL1941-0,1
 0,1 ML 1/100 100 TEST

POSITIVE CONTROL: Colon Carcinoma

VOLUME: 7 ml Ready to Use (7 ml of antibody prediluted in 0.05mol/L Tris-HCl, pH 7.6 containing

stabilizing protein and 0.015mol/L sodium azide.)

**HOST**: Mouse

CLONE: ZR347

**ANTIBODY CONCENTRATION:** 200ug/ml

SPECIES REACTIVITY: Human. Others-not known.

**EPITOPE:** Not determined

**MICROBIOLOGICAL STATE:** This product is not sterile.

**PRETREATMENT:** Staining of formalin-fixed tissue sections requires treating the tissue sections in boiling 10mM citrate buffer, pH 6.0, for 10-20 minutes followed by cooling at room temperature for 20 min.

**PRIMARY ANTIBODY INCUBATION TIME:** 30 minutes at Room Temperature

**STAINING TIPS:** If the staining is too light, use lower dilution or longer time. If the staining is

too strong, check pretreatment, use higher dilution or shorter time.

STORAGE AND STABILITY: This product contains sodium azide and is stable for 24 months when stored

at 2-8°C. Do not use after expiration date indicated on label of the product. If reagent

is not stored as recommended, performance must be validated by the user.

**TROUBLESHOOTING:** Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).



