



Microscope Cover Glass 24x50, 24x60

INTENDED USE:

For In Vitro Diagnostic Use

DESCRIPTION

Microscope Cover Glasses are used to coverslip slides by using mounting mediums. Our Cover Glasses can be used with xylene or water based mounting mediums to cover on Standard, poly-lysine or positive Charged microscope slides.

USAGE

After dropping 1-2 drops on the slide, put the cover glass on the slides. It will keep each other.

IDENTIFICATION OF THE PRODUCT AND COMPANY

Trade name: Microscope Cover Glass 24x50, 24x60

Supplier: PatoLab Biyomedikal Cihazlar ve Sarf Malzemeler Sanayi ve Tic. Ltd. Şti.

Product code: 165016, 165017

Use: Cover Glass for coverslipping slides.

HAZARDS IDENTIFICATION

It does not contain any hazard.

FIRST-AID MEASURES

Inhalation: -

Ingestion: -

Skin contact: -

Eye contact: -

FIRE FIGHTING MEASURES

Flash point: -

Ignition Temperature: -

Flammable Limit (% by vol.): -

Extinguishing Media: -

PHYSICAL AND CHEMICAL PROPERTIES

Appearance: glass

Colour: white

Odour: -

Water solubility (at 20°C): 0%;

Boiling point: -

Melt point: -

Viscosity (at 20°C): -

ph: -

Temperature at use: normal temperature;

Origin Country : Peoples Republic of China

Troubleshooting

Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).