



## PatoCool Rapid Freezing Spray

### **Intended Use :**

For In Vitro Diagnostic Use

### **Description :**

PatoCool Rapid Freezing Spray is a spray to freeze parafin tissue blocks down to - 50 C degrees on Microtome, Cryostat instruments.

### **Usage :**

Before handling be careful. Insert the tip, aim the block and push the top of the spray.

### **Identification of the Product and Company :**

**Trade name:** PatoCool Rapid Freezing Spray

**Supplier:** Abroad

**Product code:** PL-200-FRZ

**Use:** To freeze parafin tissue blocks on Microtome, Cryostat instruments.

### **Hazards Identification :**

No hazard.

### **Fire Fighting Measures :**

**Flash point:** Yes

**Ignition Temperature:** -

**Flammable :** Yes

**Extinguishing Media:** Yes

### **Physical and Chemical Properties :**

**Appearance :** Gas

**Colour :** -

**Odour :** -

### **Troubleshooting :**

Please contact Patolab Technical Support by e-mail ([patolab@patolab.com.tr](mailto:patolab@patolab.com.tr)).