

Immersion Oil

INTENDED USE : For In Vitro Diagnostic Use

<u>Catalog No</u>	<u>Volume</u>
PL-050-IMS	50 ml
PL-125-IMS	125 ml
PL-500-IMS	500 ml
PL-1000-IMS	1000 ml

DESCRIPTION

Immersion oil for optical microscope.

STORAGE & SHELF LIFE

Store at room temperature. Stable for minimum 3 years.

MICROBIOLOGICAL STATE

Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED

N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE

- 1) Place a drop of oil at the area of the slide to be examined.
- 2) Then slowly approach with the immersion lens until you touch the oil and focus through the micrometric screw.

Troubleshooting

Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).