

Paraffin Wax

INTENDED USE : For In Vitro Diagnostic Use

AÇIKLAMA :

It is produced to be used in Tissue Embedding and Tissue Processing processes.

It is white in color and in bead form.

It has a material compatible with cassettes and does not separate from the cassette.

We have 2 types of paraffin. It has melting degrees of 56-58 and 58-60.

Code	Melting Point	Package
Pwax56	56-58°C	10 Kg
Pwax58	58-60°C	10 Kg

Package Dimensions :

270mm x 270mm x 270 mm

Storage & Shelf life :

Room Temperature, 5 years

Microbiological State :

Product(s) are not sterile.

Troubleshooting :

Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).