



PatoSLim Super Mega Cassette (Slim)

Intended Use :

For In Vitro Diagnostic Use

Description :

These are lidded cassettes designed to be used for large tissue pieces in the Tissue Tracking process.

PatoSLim Super Mega Cassette can be used for large tissue pieces.

There is a compartment that allows 2 different large pieces of the same patient to be placed in the cassette.

Patient and part information can be written in the writing area on the front of the cassette.

Cassettes are resistant to high temperatures. There are 3 different Main Color options.

Code	Color	Package
PSLC100-W	White	100 pcs / bag
PSLC100-B	Blue	100 pcs / bag
PSLC100-P	Pink	100 pcs / bag

Package Dimensions :

76mm x 53mm x 7 mm

Storage Conditions :

Room Temperature

Microbiological State :

Product(s) are not sterile.

Troubleshooting :

Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).