

AQUEOUS MOUNTING MEDIA

INTENDED USE:

For In Vitro Diagnostic Use

<u>Catalog No</u>	<u>Volume</u>
PL-060-AQU	60 ml
PL-125-AQU	125 ml
PL-500-AQU	500 ml
PL-1000-AQU	1000 ml

DESCRIPTION

Aqueous Mounting Media is formulated for covering tissue sections and cell preparations immunohistochemically stained with chromogens producing alcohol soluble OR alcohol insoluble end products. Aqueous Mounting Media can be used with chromogens such as Fast Red or AEC that require an aqueous mounting medium, or those that do not, such as DAB. No heating is required prior to use. **This product is optimal for use both with and without a glass coverslip.**

WARNINGS & PRECAUTIONS

Refer to MSDS.

STORAGE & SHELF LIFE

Store at room temperature, stable for 24 months.

MICROBIOLOGICAL STATE

Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED

N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE

Apply several drops of Aqueous Mounting Media to stained specimens. Allow to dry.

Troubleshooting

Please contact PatoLab Technical Support by e-mail (patolab@patolab.com.tr).