

PatoFix Smear Spray

INTENDED USE : For In Vitro Diagnostic Use

Catalog No

PL-125-FIX

Volume

125 ml

DESCRIPTION

PatoFix is used to fix the cytology smear samples.

STORAGE & SHELF LIFE

Store at room temperature. Stable for minimum 3 years.

MICROBIOLOGICAL STATE

Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED

N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE

1) Spray 10 cm away to cytological sample on the slide.

Troubleshooting

Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).