

IN VITRO DIAGNOSTIC DATASHEET

PatoGlue Cytology Cell Sticker

INTENDED USE:

For In Vitro Diagnostic Use

Catalog No

PLC-100

Volume

100 ml

DESCRIPTION

PatoGlue is dropped into the general cytology preservative to capture the important cells in it. It captures and precipitates important cells after vortexing and centrifugation.

WARNINGS & PRECAUTIONS

Refer to MSDS.

STORAGE & SHELF LIFE

Store at room temperature, stable for 24 months.

MICROBIOLOGICAL STATE

Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED

N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE

Apply 2 drops to the general cytology preservative. After vortex and centrifugation, important cells will be captured. Apply cells to the microscope slide.

Troubleshooting

Please contact PatoLab Technical Support by e-mail (patolab@patolab.com.tr).