



Positive Charged + Microscope Slides

INTENDED USE:

For In Vitro Diagnostic Use

DESCRIPTION

Positive Charged Microscope Slide is used for immunohistochemistry procedures. It is better than Standard slides. Positive charge catches the tissue and decrease the possibility of tissue spilling.

USAGE

Get the tissue from the water bath directly. Be sure that expiry date still exist.

IDENTIFICATION OF THE PRODUCT AND COMPANY

Trade name: Positive Charged + Microscope Slide

Supplier: PatoLab Biyomedikal Cihazlar ve Sarf Malzemeler Sanayi ve Tic. Ltd. Şti.

Product code: 165013

Use: Special Microscope Slide for immunohistochemistry use.

HAZARDS IDENTIFICATION

It does not contain any hazard.

FIRST-AID MEASURES

Inhalation: -

Ingestion: -

Skin contact: -

Eye contact: -

FIRE FIGHTING MEASURES

Flash point: -

Ignition Temperature: -

Flammable Limit (% by vol.): -

Extinguishing Media: -

PHYSICAL AND CHEMICAL PROPERTIES

Appearance: glass

Colour: white

Odour: -

Water solubility (at 20°C): 0%;

Boiling point: -

Melt point: -

Viscosity (at 20°C): -

ph: -

Temperature at use: normal temperature;

Origin Country : Peoples Republic of China

Troubleshooting

Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).