



PatoCool Freezing Spray

Intended Use :

For In Vitro Diagnostic Use

Description :

PatoCool Rapid Freezing Spray is a spray to freeze parafin tissue blocks down to - 50 C degrees on Microtome, Cryostat instruments.

Usage :

Before handling be careful. Insert the tip, aim the block and push the top of the spray.

Identification of the Product and Company :

Trade name: PatoCool Rapid Freezing Spray

Supplier: Abroad

Product code: PL-200-FRZ

Use: To freeze parafin tissue blocks on Microtome, Cryostat instruments.

Hazards Identification :

No hazard.

Fire Fighting Measures :

Flash point: Yes

Ignition Temperature: -

Flammable : Yes

Extinguishing Media: Yes

Physical and Chemical Properties :

Appearance : Gas

Colour : -

Odour : -

Troubleshooting :

Please contact Patolab Technical Support by e-mail (patolab@patolab.com.tr).